

# **Structure Determination Tutorial**

Michael Baran  
G.V.T. Swapna  
Jim Aramini  
Roberto Tejero  
Janet Huang  
Nadeem Riaz

6/26/03

7/28/03

## Contents

I. Process Backbone Experiments using standard AutoProc parameters

II. Peak Picking

II.A Identify NH<sub>2</sub> Residues in HSQC

II.B Identify folded peaks from large sw HSQC

III.B Peak Pick Triple Resonance Experiments

III. Backbone Assignments

IV. Side Chain Assignments

IV.A Aliphatic Sidechain assignments

V.B Assign NH<sub>2</sub>'s

IV.C Aromatic Side Chain Assignments

IV.D Run missing\_shifts.pl

V. Structure Determination

V.A. AutoStructure

V.B Run Talos for dihedral angle constraints

V.C Analyze HNHA experiments to determine J values

V.D Xplor Refinement

## **I.A Process Backbone Experiments using standard AutoProc parameters**

Process each experiments xy and xz plane to set correct phasing.

Use proper zero filling. CABM data collection usually use 1024 x 512 x 256 matrix. For C13-edited NOESY use 1024x512x512

Use | nmrPipe -fn EXT -left -sw flag for 3d processing where appropriate. This command will exclude half of spectra where no peaks appear and thus reduce processing time.

## **II. Peak Picking (Triple Resonance Experiments)**

### **II.A Identify NH2 Residues in HSQC**

1. Open HSQC experiment.
  1. Display what nuclei is on each axis
    - i. Sparky cmd: xa
  2. Rotate cube so x1 is H and x2 is N
    - i. Sparky cmd: xx and xr
  3. Set proper contour levels
    - i. Sparky cmd: ct
  4. Adjust aspect ratio
    - i. Sparky cmd: vt
2. Overlay NH2-HSQC onto NH-HSQC
  1. Sparky command - ol
3. Pick all HSQC peaks
  1. Sparky cmd: F11 and draw box around all peaks to select
4. Remove NH2 peaks from HSQC
  1. Select peaks with F1 and delete from list
5. Run examine\_expected\_peaks.pl to be sure all peaks are accounted for OR manually determine how many peaks expected
  1. Num peaks = num or residues - # of PRO - 1 terminal NH?

### **II.B Identify folded peaks from large sw HSQC**

### **III.B Peak Pick Triple Resonance Experiments**

1. Peak pick HNCO, HNCACB, HNcoCACB and any other triple resonance experiments against HSQC.
2. Open triple resonance experiment as described in II.A.1
3. Synchronize spectra to HSQC
  - a. Sparky cmd: yt
  - b. Synchronize each nuclei individually
4. Restrictive peak pick spectra
  - a. Sparky cmd: kr
  - b. use default tolerances
  - c. adjust contour levels (ct) as needed until roughly expected number of peaks are picked
    - i. CO (intra) experiments – expect approx same number of peaks as residues
    - ii. no CO (inter) experiments – expect 2X as many peaks as residues
5. Convert peak lists to autoassign format and perform peak list validations
  - a. run extract\_columns.pl – to extract chemical shift columns and data height
  - b. run create\_peak\_list.pl – to convert extract\_columns output to autoassign format lists
6. Run examine\_expected\_peaks on each spectra to evaluate if number of peaks picked is reasonable
7. Run calculate\_registration.pl – large registration errors can be used to detect flawed peak lists.
  - a. register each spectra to HSQC to register N H dimensions
  - b. register C dimension by registering spectra to HNcoCACB
8. Run examine\_spin\_systems.pl report. ESS report can be used to detect incomplete or overlapped spin systems.

### **III. Backbone Assignments**

1. Convert peak lists to AutoAssign format as described above in step II.B.5.
2. Register peak lists using calculate\_registration.pl as described in step III.B.7
  - a. Generate strategy specific peak list (autoassign runs better with HNCA and HNcoCA data. However, if these experiments are absent they can be “created” from HNCACB and HNcoCACB experiments
    - i. Ex. Gen HNCA from HNCACB
    - ii. Ex2. Generate HNcoCA from HNcoCACB
    - iii. Use extract\_by\_phase and extract\_by\_filter.pl to produce these “pseudo-peak lists”
  - b. Use registration value to shift peaks lists using apply\_shifting.pl
  - c. Calculate AutoAssign tolerances
    - i. Take largest Standard deviation and apply by 3 or 4 for each dimension
    - ii. Example of typical tolerances are: N: 0.01 H: 0.11 CA: 0.21

CB: 0.21 CO: 0.4 HA: 0.02

3. Write AutoAssign Control file
4. Write AutoAssign Override file
  - a. overrides from HCCH-TOCSY
  - b. Gly phased HA experiment overrides
  - c. Selective perdeuteration overrides
5. Run AutoAssign to make assignments
  - a. Save cmap files/images of assignments
  - b. Save bmrB file of assignments
6. Run validate\_assignments.pl to evaluate problematic assignments
7. Paste assignments back into spectra for manual validation
  - a. BmrB2SparkyPks.pl - script to convert autoassign output bmrB file to a peak list which can be opened in Sparky.
  - b. Open sparky HSQC
  - c. Read Assigned Sparky List from BmrB2SparkyPks.pl into each peak list
    - i. Sparky cmd: rp
  - d. Synchronize Peak lists and verify assignments.
  - e. Edit peak lists and re-run autoassign until satisfied.

#### **IV. Side Chain Assignments**

##### **IV.A Aliphatic Sidechain assignments**

1. Process HCCH-TOCSY and hCCH-TOCSY
2. Use these spectra with triple resonance backbone experiments and average carbon and hydrogen side shift values to tentatively assign all spin system Carbons and Hydrogens.
3. Process HCCH-COSY
4. Make HH and CH projection of HCCH-COSY
  - a. Sparky cmd: (XH) ucsfdata -s1 512 -r -o [out] [in]
  - b. Sparky cmd: (HH) ucsfdata -s1 256 -r -o [out] [in]
  - c.
5. Use projections to restrictive peak HCCH-COSY
6. Run 1 cycle of AutoStructure to assign HCCH-COSY
7. Paste AutoStructure assignments into Sparky HCCH-COSY
  - a. Sparky cmd: rp
8. Use AutoStructure QM file to help find un-assigned residues
9. Manually confirm side chain assignments using both HCCH-TOCSY ,hCCH-TOCSY and HNCACB.
  - a. Use expected H chemical shift chart as well as expected peak pattern knowledge as a guide.
  - b. Frequently view assignments standard deviations using the sparky rl command. Keep all standard deviations below 0.1.
10. Validate side chain assignments using validate\_assignments.pl
  - a. Save rl results from sparky

- b. Use createBmrb.pl to generate bmrB file used as input for validate\_resonance\_Assignments.pl

#### **IV.B Assign NH<sub>2</sub>'s**

I would suggest assigning NH's after running an initial AutoStructure run. This way the assignments from AutoStructure can be pasted into the NOESY to locate the given residue strip.

Use NOESY's to stereo-specifically assign NH<sub>2</sub>'s. based on peak intensities and expected peak positions.

#### **IV.C Aromatic Side Chain Assignments**

Again I would suggest that aromatic assignments are not made until an initial AutoStructure run has been completed. This way assignments can be pasted into the C13-edited NOESY and N15-edited NOESY.

Process Aromatic C13-edited HSQC and aromatic NOESY.

Use these experiments along with aliphatic NOESY's to assign aromatics.

(If using Yamasaki RD experiments good luck!)

#### **IV.D Run missing\_shifts.pl**

Run AutoPeak missing\_shifts.pl to evaluate assignment completeness.

### **V. Structure Determination**

#### **V.A. AutoStructure**

NOTE: use -j flag to include jcoupling data

HYPER and Talos may conflict w/ each other, use caution when running both

1. Process N15-edited NOESY and C13 edited NOESY
  - a. Make sure to use LP!! Process C13-NOESY w/ 1024x512x512 matrix
2. Peak pick noesy's against respective HSQC as well as HH projection.
3. Generate AutoStructure Input
4. Run auto-structure
5. analyze AutoStructure QM file for missed residues in peak lists
6. Re-run auto structure until satisfied with quality scores.
  - a. M score < 0.1 N15 NOESY
  - b. M score < 0.4 C13 NOESY
7. Then run multiple cycles of AS analyzing the unassigned file.

8. Remove spurious peaks which are being unassigned and re-run AS until ~under 10 assignments are being removed per cycle.
9. Calculate I,L scores

### **V.B Run Talos for dihedral angle constraints**

1. run Tab4Talos.pl script to convert bmrB assignments into Talos input
2. run talos
3. run Talos2Dyana.pl to convert talos output to Dyana format for use in AutoStructure

### **V.C Analyze HNHA experiments to determine J values**

NOTE: Jvalues 6-8 sheet

2-4 values

1. peak pick HNHA
2. In order to calculate J values need to separate diagonal and off-diagonal peaks
  - a. **OFF-DIAG peaks**
  - b. BmrB2SparkyPks.pl HA N H
  - c. Paste rp
  - d. BmrB2SparkyPks.pl HA2 N H
  - e. Paste rp
  - f. BmrB2SparkyPks.pl HA3 N H
  - g. Paste rp
  - h. **DIAGONAL PEAKS**
  - i. BmrB2SparkyPks.pl H N H
  - j. Paste rp
3. Create SpreadSheet in Excel
4. columns
  - a. Diag Assignment
  - b. Diag peak intensity
  - c. Off- Diag peak assignment
  - d. Off-Diag peak intensity
  - e.  $1.1 * I(x) / I(d)$
  - f.  $\sqrt{1.1 * I(x) / I(d)}$
  - g.  $\text{atan}(\sqrt{1.1 * I(x) / I(d)})$
  - h.  $\text{atan}(\sqrt{1.1 * I(x) / I(d)}) / 2 * \text{pie} * 0.0125$ 
    - i. excel pie = 3.14 cladiograph pie = 180
5. Convert to bmrB format manually

assign	diag	off-diag ass	off-diag I
I3H-N-H	405877 I3HA-N-H	-146432	-0.396857176 0.629966012 0.562162412 7.161304616
I4H-N-H	164581 I4HA-N-H	-125338	-0.837713952 0.915267148 0.741186245 9.441862988
S5H-N-H	823571 S5HA-N-H	-161818	-0.2161317 0.464899666 0.435175174 5.543632785

I6H-N-H	564057 I6HA-N-H	-301423	-0.58782233	0.766695722	0.654101023	8.332497114
S7H-N-H	351641 S7HA-N-H	-211607	-0.661946986	0.813601245	0.682979536	8.700376249
T9H-N-H	428115 T9HA-N-H	-273439	-0.702575009	0.838197476	0.697602059	8.886650437
N11H-N-H	546969 N11HA-N-H	-346288	-0.696413874	0.834514154	0.69543475	8.859041395
N13H-N-H	347744 N13HA-N-H	-151580	-0.479484908	0.692448488	0.605639835	7.715157136
T14H-N-H	369613 T14HA-N-H	-213263	-0.634688986	0.796673701	0.672709417	8.569546719
M15H-N-H	246989 M15HA-N-H	-191134	-0.851241958	0.922627746	0.745176895	9.492699302
K16H-N-H	315093 K16HA-N-H	-175804	-0.613737531	0.783414023	0.664545403	8.465546532
I17H-N-H	168571 I17HA-N-H	-103481	-0.67525909	0.821741498	0.687858078	8.762523283
T18H-N-H	314057 T18HA-N-H	-152424	-0.533872514	0.730665802	0.631011961	8.038368926
L19H-N-H	167629 L19HA-N-H	-107724	-0.706896778	0.840771537	0.699112008	8.905885455
S20H-N-H	454340 S20HA-N-H	-55344	-0.133993045	0.366050604	0.350901634	4.470084507

J BMRB format for AutoStrucutre

1 3JHNHA	3 ILE	H	3 ILE	HA	7.161305	1.5
2 3JHNHA	4 ILE	H	4 ILE	HA	9.441863	1.5
3 3JHNHA	5 SER	H	5 SER	HA	5.543633	1.5
4 3JHNHA	6 ILE	H	6 ILE	HA	8.332497	1.5
5 3JHNHA	7 SER	H	7 SER	HA	8.700376	1.5
6 3JHNHA	9 THR	H	9 THR	HA	8.88665	1.5
7 3JHNHA	11 ASN	H	11 ASN	HA	8.859041	1.5
8 3JHNHA	13 ASN	H	13 ASN	HA	7.715157	1.5
9 3JHNHA	14 THR	H	14 THR	HA	8.569547	1.5
10 3JHNHA	15 MET	H	15 MET	HA	9.492699	1.5
11 3JHNHA	16 LYS	H	16 LYS	HA	8.465547	1.5
12 3JHNHA	17 ILE	H	17 ILE	HA	8.762523	1.5
13 3JHNHA	18 THR	H	18 THR	HA	8.038369	1.5
14 3JHNHA	19 LEU	H	19 LEU	HA	8.905885	1.5
15 3JHNHA	20 SER	H	20 SER	HA	4.470085	1.5
16 3JHNHA	21 GLU	H	21 GLU	HA	8.226777	1.5

## V.D Xplor Refinement

1. Use following script:
  - a. ShifttoXplor <assignment.bmr> [output]
2. Use Xplor constraints from AutoStructure
  - a. Read xplor constraints into PdbStat and rename nomenclature conflicts
    - i. Rea pdb XX.pdb
    - ii. Rea cons xplor XX.xpl
      1. find bad Me-Me contacts
      2. change all CZ to HE#
      3. change all CG to HD#
3. Convert DYANA aco file to Xplor w/ PdbStat (may want to remove HYPER constraints and keep only Talos constraints before conversion to Xplor format)
  - a. Rea aco xplor XX\_dihe.tbl
4. Convert Hbonds to Xplor format w/ PdbStat



- a. rea upl dyana hbond.upl
- b. rea lol dyana hbond.lol
- c. write cons xplor XX\_hbond.tbl
5. generate topology file
  - a. /farm/software/xplor-nih-2 <generate\_psf
6. make pdb for Xplor
  - a. rea coo pdb XX.pdb
  - b. write coor pdb XX\_0.pdb
7. ./RefineProc XX 10 10 PBS 0
8. Make PDB file of 10 best structures
  - a. Concatenate all 100 pdb files
  - b. Agrupa zr18\_sa\*.pdb
  - c. Pdbstat
  - d. Rea coors all
  - e. Classify
  - f. Del 11 100
  - g. Order 0.9
  - h. Rmsd best backbone
  - i. Write coor pdb best.pdb all

NOTES on getting Xplor to work:

1. Be sure you are added as a user to PBS with proper permissions
2. be sure paths are set correctly for DQS or PBS in RefineProc script
3. set the following paths
  - a. export TOPPAR=/farm/software/xplor-nih-2.0.6/toppar
  - b. export PATH=/farm/software/bin:\${PATH}
  - c. export C13SHIFTS=/farm/software/xplor-nih-2.0.6/databases/c13shifts
  - d. export GAUSSIANS=/farm/software/xplor-nih-2.0.6/databases/torsions\_Gaussians

## **Ribbons**

### **GENERAL EXECUTION**

Creating a Protein structure info

1. run ribbons-model <pdb> <model name> Produces ss structure info
2. Run atom-model <pdb> Produces atom sphere and bond info
3. To create disulfide bond info, go to File->Export->SG Bonds
4. Create ribbon to side chain bond info File->Export->CG Bonds
5. Polish CG bonds so only have active component parts

Coloring Secondary Structures

Add another tag to

Important programs

pdb-range-ss: allows for coloring certain ss different colors than default  
pdb-sele-pdb: allow to cut certain residues outa pdb file

### **MODEL CREATION**

Steps:

1. get PDB
2. Created base model structure (ribbons-model <pdb file>)
3. Create atom model (atom-model <pdb file>)
4. Correct secondary structure to match pdb file
5. If certain important functional sites, create pdb file w/only those residues
  
7. run atom-model on aromatic pdb file
8. run pdb-res-sph with color.1 color file to get color coded sphere
9. run sph-bond on file created in 8 to get color coded cylinders (coded by residue) \* Note 8 & 9 need to be rerun to change colors of cylinders \*
10. Done!
11. correct \_cg\_bonds.cyl file

For texts, see RBFA for example .text file

### **RIBBONS CODES**

Ribbons Coloring Codes:

#	Color
1	Red
2	Green
3	Yellow
4	Blue
5	Magenta
6	Cyan
7	White
8	Orange
9	Grey
10	Lapis
11	Silver
12	Lavendar
13	Brown
14	Gold
15	Rose
16	Lightgreen

Ribbons Secondary Structure Codes

Letter      Structure

S Sheet  
A Sheet w/Arrow  
H Right handed alpha helices  
L Left handed alpha helices  
3 3-10 helices  
c coil

### **RIBBONS FILE INFO**

.ribbons file contains names of polypeptide ss files and optional color file

.ss file contain polypeptide ss info

First line is arbitrary title line

Second line contains coloring keys, 3 or 4 mandatory keys, rest are user definable

.colors file contain color info

See Ribbons-CODES for color code info

.sph contains info on set of atoms (i.e. geometrical representation)

.cyl contains info on set of bonds (i.e. geometrical representation)

.atoms contains list of .sph files (so you can specificall alter one set of atom s)

.bonds contains list of .cyl files

\_sg.cyl disulfide bond info

\_cb.cyl create bonds from side chains to ribbons

the name field (last field) contains residue info, use that to determine bonds from ribbons to keep and which to dum

p

### **PYMOL**

Very simple to use.

1. Be sure to add secondary structure information to header of PDB file
  - a. /farm/software/structure\_validation/get\_sse\_info.pl
  - b. cut and paste results of this script into header of pdb
  - c. open pdb file with pymol and have fun.

## **GRASP**

The following is how to open your molecule, create a charged surface and create 2 views front and side for publication using grasp. Use the right mouse button to bring up menu as well as select

Calculate charge of molecule:

Read -> Radius/charge File (+assign) -> full.crg

Build Surface:

Build -> Molec Surface -> all atoms

Calculate potential:

Calculate -> New Potential Map -> all atoms

Calculate -> Pot. Via Map and Surface -> Atoms -> all atoms -> all surface

Then adjust values on top bar to increase/decrease redness

To Duplicate images (for front and side picture)

Display -> stereo/split on -> dials both

To bring up backbone wire view:

You must build the backbone wire view. (sorry no command here) Then to see through the potential surface

Panel -> objects -> surface -> transparent